

Documentation of data processing

This file contains the documentation of the scripts in the `analysis/scripts` folder. The documentation has been written in literate markdown. To produce the scripts from the source file (`process.praat.md`), use the `lmt` package (written in Go by Dave MacFarlane, at <https://github.com/driusan/lmt>).

Prepare files for force alignment in SPPAS

`alignment-input.praat`

```
<<<get audio>>>

<<<concatenate recoverably>>>

<<<write sppas>>>
```

The following chunk asks the user for the name of the project directory and the participant ID. Then it reads the audio files from the `audio` directory (which contains the audio files exported from AAA). The directory `alignment` is created as well.

“get audio”

```
form Generate input for force alignment with SPPAS
  word project pilot
  word speaker sc01
endform

directory_speaker$ = "../data/derived/ultrasound/'speaker'"
directory_audio$ = "'directory_speaker$'/audio"
createDirectory ("'directory_speaker$'/alignment")
directory_alignment$ = "'directory_speaker$'/alignment"
writeFile: "'directory_alignment$'/speaker$.txt", ""

Create Strings as file list: "filelist", "'directory_audio$'/*.wav"
files = Get number of strings

for file from 1 to files
  select Strings filelist
  file$ = Get string: file
  Read from file: "'directory_audio$'/'file'"
  sound = selected("Sound")
  sound$ = file$ - ".wav"
endfor
```

The following select all objects except the file list and concatenates the sound objects. SPASS needs a tier named `Orthography`, so we create one. Then, the script loops through the intervals in the

TextGrid which correspond to the names of the files, it writes the prompt in the interval and in the text file for IPU detection.

“concatenate recoverably”

```
select all
minusObject: "Strings filelist"
Concatenate recoverably

selectObject: "TextGrid chain"
Duplicate tier: 1, 1, "Orthography"

intervals = Get number of intervals: 1

for interval from 1 to intervals
  start = Get start point: 1, interval
  end = Get end point: 1, interval
  filename$ = Get label of interval: 1, interval

  Read Strings from raw text file: "'directory_audio$/'filename$.txt"
  prompt$ = Get string: 1
  selectObject: "TextGrid chain"
  Set interval text: 1, interval, "'prompt$"
  appendFileLine: "'directory_alignment$/'speaker$.txt", "'prompt$"
endfor
```

Finally, we can save the concatenated sound file and the TextGrid with the file names. The latter will be used in the script `search-area.praat` to separate the concatenated TextGrid.

“write sppas”

```
selectObject: "Sound chain"
Save as WAV file: "'directory_alignment$/'speaker$.wav"

selectObject: "TextGrid chain"
Copy: "filenames"
Remove tier: 1
Save as text file: "'directory_alignment$/'speaker$-filenames.TextGrid"
```

Extract the search area for spline batch processing and kinematics in AAA

`search-area.praat`

```
<<<get alignment>>>

<<<set search>>>
```

```
<<<extract search>>>
```

The user is prompt to indicate the project name, the participant ID and the language. Depending on the language selected, the appropriate speech segments are stored for subsequent extraction of the search area. Then the script read the TextGrid file with the force alingment (`ID-align.TextGrid`). The number of intervals of the TextGrid file is saved in `intervals` and two new tiers are created (ultrasound and kinematics)

“get alignment”

```
form Select folder with TextGrid
  word project pilot
  word speaker sc01
  comment Supported languages: it, pl
  word language it
endform

if language$ == "it"
  label_lang$ = "k"
  label_2_lang$ = "dico"
elif language$ == "pl"
  label_lang$ = "j"
  label_2_lang$ = "mowie"
else
  exit "The language you selected is not valid"
endif

directory_audio$ = "../data/derived/ultrasound/'speaker$'/audio"
directory_alignment$ = "../data/derived/ultrasound/
... 'speaker$'/alignment"

palign = Read from file: "'directory_alignment$'/'speaker$'-palign.TextGrid"
palign.o = Read from file: "'directory_alignment$'/'speaker$'-palign.TextGrid"
selectObject: palign

intervals = Get number of intervals: 1
```

Now we can create intervals containg the search area for ultrasound and kinematics which will be used in AAA for spline batch processing and to find consonantal gestures moments. Then, `[ID]-search.TextGrid` is saved in the `alignment` folder.

“set search”

```
Insert interval tier: 4, "ultrasound"
Insert interval tier: 5, "kinematics"
Insert interval tier: 6, "vowel"

for interval to intervals
```

```

label$ = Get label of interval: 1, interval
if label$ == label_lang$
    start_ultrasound = Get start time of interval: 1, interval
    interval_2 = Get interval at time: 2, start_ultrasound
    label_2$ = Get label of interval: 2, interval_2
    if label_2$ == label_2_lang$
        end_ultrasound = Get end time of interval: 1, interval + 7
        Insert boundary: 4, start_ultrasound
        Insert boundary: 4, end_ultrasound
        ultrasound = Get interval at time: 4, start_ultrasound
        Set interval text: 4, ultrasound, "ultrasound"

        start_kinematics = Get start time of interval: 1, interval + 3
        end_kinematics = Get end time of interval: 1, interval + 5
        start_vowel = Get start time of interval: 1, interval + 3
        end_vowel = Get end time of interval: 1, interval + 3

        if label_2$ == "dico"
            start_kinematics = ((end_vowel - start_kinematics) / 2) +
                ...start_kinematics
            start_vowel_2 = Get start time of interval: 1, interval + 5
            end_kinematics = ((end_kinematics - start_vowel_2) / 2) +
                ...start_vowel_2
        endif

        Insert boundary: 5, start_kinematics
        Insert boundary: 5, end_kinematics
        kinematics = Get interval at time: 5, start_kinematics
        Set interval text: 5, kinematics, "kinematics"

        vowel$ = Get label of interval: 1, interval + 3
        Insert boundary: 6, start_vowel
        Insert boundary: 6, end_vowel
        vowel_interval = Get interval at time: 6, start_vowel
        Set interval text: 6, vowel_interval, vowel$
    endif
endif
endfor

Remove tier: 1
Remove tier: 1
Remove tier: 1

Save as text file: "'directory_alignment$'/'speaker$'-search.TextGrid"

```

Then, the script saves each search area to separate TextGrids in the `audio` folder. The file names are extracted from `[ID]-filenames.TextGrid`.

“extract search”

```
filenames = Read from file: "'directory_alignment$'/'speaker$'-filenames.TextGrid"

selectObject: palign
plusObject: filenames

Merge

filenames_tier = 4

intervals = Get number of intervals: filenames_tier

for interval from 1 to intervals
  selectObject: "TextGrid merged"
  start = Get start point: filenames_tier, interval
  end = Get end point: filenames_tier, interval
  filename$ = Get label of interval: filenames_tier, interval

  Extract part: start, end, "no"

  Remove tier: filenames_tier
  Write to text file: "'directory_audio$'/'filename$'.TextGrid"
  Remove

  selectObject: palign.o
  Extract part: start, end, "no"
  Write to text file: "'directory_audio$'/'filename$'-palign.TextGrid"
  Remove
endfor
```

Synchronise EGG data with AAA audio data

The following chunk calls the header of the script, which is defined at the end of the documentation, and the main function.

sync-egg.praat

```
<<<sync header>>>

<<<sync function>>>
```

The script works by selecting all the files in the Object window after loading files.

“sync function”

```
<<<check objects>>>

<<<read files>>>

<<<sync>>>
```

Before running, the script checks if the objects list is empty. If not, the script exits and prompts the user to clean the objects list.

“check objects”

```
select all
number_selected = numberOfSelected ()
if number_selected > 0
    exitScript: "Please, remove the objects in the Objects window. For this
    ... script to work, the Objects list must be empty."
endif
```

The form asks for the project name and the participant ID. A boolean is stored as well for enabling the debug mode. In the debug mode, all intermediate files produced by the script are kept in the Objects window. They are deleted otherwise.

“read files”

```
form Synchronise EGG data
    word project pilot
    word speaker sc01
    boolean debug_mode
endform
```

The file lists of the EGG and ultrasound .wav files are saved in `filelist_egg` and `filelist_us`. The number of files in the EGG folder is saved in `files`.

“read files”+=

```
egg_directory$ = "../data/raw/egg"
us_directory$ = "../data/derived/ultrasound"
out_directory$ = "../data/derived/egg"
createDirectory ("out_directory$/'speaker$'")

Create Strings as file list: "filelist_egg", "egg_directory$/'speaker$'/*.wav"
files = Get number of strings
Create Strings as file list: "filelist_us", "us_directory$/'speaker$'/audio/*.wav"
```

For every file listed in `filelist_egg`, it reads the file.

“read files”+=

```
for file from 1 to files
  select Strings filelist_egg
  file$ = Get string: file
  Read from file: "'egg_directory$'/'speaker$'/'file$'"
endfor
```

Every object is then selected, minus the two file lists. The Sounds are concatenated to **Sound chain**.

“read files”+=

```
select all
minusObject: "Strings filelist_egg"
minusObject: "Strings filelist_us"
Concatenate
```

While **Sound chain** is selected, the script inverts the signal (both the audio and the EGG signal are inverted during acquisition with the Laryngograph), and extracts all channels (the object is a stereo sound: channel 1 is the audio, channel 2 is the EGG signal). For cross-correlation to work, the two sound files must have the same sampling frequency. AAA records at a frequency of 22050 Hz. To ensure that the EGG audio is at the same sampling frequency, resampling is performed. **Sound chain_ch1_22050** is created from **Sound chain_ch1**.

“read files”+=

```
Multiply: -1
Extract all channels

selectObject: "Sound chain_ch1"
Resample: 22050, 50
```

The extraction of each stimulus from the concatenated sound is achieved through the EGG signal. The function **To TextGrid (silences)** efficiently recognises the voiced stretches of the audio which roughly corresponds to the spoken stimuli. (*Warning*: this assumes that the EGG files don't contain spurious material) The minimum duration for silence is set to 1 second to avoid voiceless segments being annotated as silence. The output is **TextGrid chain_ch2**. The number of intervals in the TextGrid is saved in the variable **intervals**.

“sync”

```
selectObject: "Sound chain_ch2"

To TextGrid (silences): 100, 0, -25, 1, 0.1, "silence", "speech"
intervals = Get number of intervals: 1
Insert interval tier: 2, "new"
```

For each interval in the TextGrid `chain_ch2` which is labelled `speech`, the start and end time of the interval are moved by -1.5 and 1 second respectively. This ensures that there is enough audio before and after the stimulus for cross-correlation. The original left and right boundaries are removed. The result is that the interval label is changed to `speechsilence`.

“sync”+=

```
for interval from 1 to intervals
  label$ = Get label of interval: 1, interval
  if label$ == "speech"
    start = Get starting point: 1, interval
    end = Get end point: 1, interval
    Insert boundary: 2, start - 1.5
    Insert boundary: 2, end + 1
    new_interval = Get interval at time: 2, start
    Set interval text: 2, new_interval, "speech"
  endif
endfor

Remove tier: 1
```

We can now get the number of intervals of the updated TextGrid and set the counter `index` to 1. The counter is used to read the ultrasound audio files, and in the names of the output files.

“sync”+=

```
intervals = Get number of intervals: 1

index = 1
```

For every interval in the TextGrid it is checked if the label is `speechsilence`. The intervals with this label correspond to the individual stimuli in the concatenated EGG sound files. If the label is `speechsilence`, the script gets the start and end time of that interval.

“sync”+=

```
for interval from 1 to intervals
  label$ = Get label of interval: 1, interval
  if label$ == "speech"
    start = Get starting point: 1, interval
    end = Get end point: 1, interval
```

Then the resampled `Sound chain_ch1_22050` is selected and the portion from `start` to `end` is extracted. This portion corresponds to the TextGrid interval and, thus, to the stimulus. The sound is named `Sound chain_ch1_22050_part`.

“sync”+=


```
selectObject: "Sound chain_ch1_22050"
Extract part: start, end, "rectangular", 1, "no"
```

The counter `index` is now used to read the audio file from the ultrasound directory. Since the order of the stimuli is the same in both the EGG and ultrasound files, a counter that increases for every interval with the `speechsilence` label is sufficient. The name of the file is saved after reading and the file remains selected.

“sync”+=

```
selectObject: "Strings filelist_us"
file_us$ = Get string: index
Read from file: "'us_directory$'/'speaker$'/audio/'file_us$'"
file_us_name$ = selected$ ("Sound")
```

The extracted portion from the EGG audio channel is added to the selection. The cross-correlation between the EGG and ultrasound audio is performed. The time of maximum amplitude in the generated cross-correlated sound corresponds to the off-set between the two files.

“sync”+=

```
plusObject: "Sound chain_ch1_22050_part"

crosscorrelated = Cross-correlate: "peak 0.99", "zero"
offset = Get time of maximum: 0, 0, "Sinc70"
```

The concatenated stereo sound (or the recombined stereo if the `invert egg signal` option is active) is selected and a portion is extracted. The portion starting point corresponds to the starting point of the TextGrid interval minus the off-set obtained from the correlation. If the offset is positive (when the audio is longer than the EGG audio), silence is added at the beginning of the EGG sound. If the offset is negative (the EGG sound is longer than the audio), the extra part is deleted from the beginning of the EGG sound to match the beginning of the audio. The end point is the same as the one of the interval. (The endpoint does not matter, since timing is calculated from the beginning of the file.) The sound is finally saved in the `sync` folder.

“sync”+=

```
selectObject: "Sound chain"

start = start - offset
Extract part: start, end, "rectangular", 1, "no"
Save as WAV file: "'out_directory$'/'speaker$'/'file_us_name$'.wav"
```

If the debugging mode is off, all the intermediate files are removed. Otherwise they are kept for inspection. The index is increased by one and the TextGrid is selected for the next cycle of the for loop.

“sync”+=

```
if debug_mode == 0
  removeObject: "Sound chain_ch1_22050_part", "Sound " + file_us_name$,
  ...crosscorrelated, "Sound chain_part"
endif

index += 1
select TextGrid chain_ch2
endif
endfor
```

Extract VUV intervals

This script calculates the voiced and voiceless portions (VUV) in the synchronised EGG files based on the EGG signal.

extract-vuv.praat

```
<<<get synced egg>>>

<<<vuv>>>
```

We first read ask for the project name and the speaker ID.

“get synced egg”

```
form Extract vuv
  word project pilot
  word speaker sc01
  boolean debug_mode
endform

directory$ = "../data/derived/egg"

Create Strings as file list: "filelist", "'directory$'/'speaker$'/*.wav"
files = Get number of strings
```

Now, for each file in derived/egg, we can calculate the boundaries of the voiced and voiceless intervals in the file and save them to a TextGrid file.

“vuv”

```
for file from 1 to files
  selectObject: "Strings filelist"
  file$ = Get string: file
```

```

Read from file: "'directory$'/'speaker$'/'file$'"
filename$ = selected$("Sound")

<<<to vuv>>>

<<<save vuv>>>
endfor

removeObject: "Strings filelist"

```

To calculate voiced and voiceless intervals, we can exploit the already available function `To TextGrid (vuv)`. The channel containing the EGG signal (channel 2) is extracted, a `PointProcess` object is created from the signal, and finally the `vuv` function is applied.

“to vuv”

Extract one channel: 2

To `PointProcess` (periodic, cc): 75, 600

To `TextGrid` (vuv): 0.02, 0.001

The resulting `TextGrid` is saved in the same synced EGG files folder.

“save vuv”

```

Write to text file: "'directory$'/'speaker$'/'filename$'-vuv.TextGrid"

if debug_mode == 0
  removeObject: "Sound " + filename$, "Sound " + filename$ + "_ch2",
  ... "PointProcess " + filename$ + "_ch2", "TextGrid " + filename$ + "_ch2"
endif

```

DEGG tracing

degg-tracing.praat

```

<<<smoothing>>>

<<<get files list>>>

<<<file loop>>>

```

First we get the file list and we start looping through the files.

“get files list”

```
form dEGG tracing
  word project pilot
  word speaker sc01
  comment Specify the lower and upper frequency (in Hz) for filtering:
  real lower 40
  real upper 10000
  comment Specify the smooth width "m" (the number of points):
  real smooth_width 11
endform

directory$ = "../data/derived/egg/'speaker$'"
directory_textgrid$ = "../data/derived/ultrasound/'speaker$'/audio"

result_file$ = "../results/'speaker$'-degg-tracing.csv"
header$ = "speaker,file,date,word,time,rel.time,proportion,maximum,minimum"
writeFileLine: "'result_file$'", "'header$'"

Create Strings as file list: "filelist", "'directory$'/*.wav"
files = Get number of strings
```

For each file, extract both channels. Read from the corresponding TextGrid in /data/derived/ultrasound/ID/audio and get the starting and end point of the vowel interval. Now, we can extract the same interval from channel 2 of the EGG file. Rename the extracted part as egg, and execute the main function, which extracts the dEGG trace.

“file loop”

```
for file to files
  selectObject: "Strings filelist"
  file$ = Get string: file
  filename$ = file$ - ".wav"

  Read Strings from raw text file: "'directory_textgrid$'/'filename$'.txt"
  prompt$ = Get string: 1
  stimulus$ = extractWord$(prompt$, " ")
  date$ = Get string: 2

  Read separate channels from sound file: "'directory$'/'file$'"

  Read from file: "'directory_textgrid$'/'filename$'.TextGrid"

  start = Get starting point: 3, 2
  end = Get end point: 3, 2

  selectObject: "Sound 'filename$'_ch2"
  Extract part: start, end, "rectangular", 1, "yes"
  Rename: "egg"
```

```
<<<main function>>>
endfor
```

“main function”

```
<<<degg>>>

<<<degg loop>>>
```

The EGG signal file `egg` is selected. Filter EGG signal (`egg_band`) and smooth it with moving average (renamed to `egg_smooth`). Create PointProcess (peaks) for EGG (`PointProcess egg_smooth`). Calculate dEGG and smooth it (?) (`degg_smooth`). Create PointProcess (peaks) of dEGG (`PointProcess degg_smooth`).

“degg”

```
Filter (pass Hann band): lower, upper, 100
@smoothing: smooth_width
Rename: "egg_smooth"
To PointProcess (periodic, peaks): 75, 600, "yes", "no"

selectObject: "Sound egg_smooth"
Copy: "degg"
Formula: "self [col + 1] - self [col]"
@smoothing: smooth_width
Rename: "degg_smooth"
To PointProcess (periodic, peaks): 75, 600, "yes", "no"
```

Loop through the EGG points and get minimum between the first two points. The loop needs to go to the number of points minus 2 since we are selecting three points in each cycle of the loop. This will need to be fixed if we want all cycles to be included. Get dEGG maximum on the right of EGG minimum and get minimum of dEGG between current maximum and the next. Normalise max and min to unity. This is gonna be the y axis. The x axis needs to be time aligned: can choose between several (use minimum in EGG as arbitrary epoch, or midway between minima, or what else?) Go to the second and third point and repeat.

ATTENTION! You don't need to normalise to unity. You need to get proportion $(\text{period} - \text{value})/\text{period}$.

Trying `egg_minimum_2` instead of `degg_maximum_2` for cases when there is no `degg_maximum_2`.

“degg loop”

```
selectObject: "PointProcess egg_smooth"
egg_points = Get number of points
mean_period = Get mean period: 0, 0, 0.0001, 0.02, 1.3
```

```

for point to egg_points - 2
  selectObject: "PointProcess egg_smooth"
  point_1 = Get time from index: point
  point_2 = Get time from index: point + 1
  point_3 = Get time from index: point + 2
  selectObject: "Sound egg_smooth"
  egg_minimum_1 = Get time of minimum: point_1, point_2, "Sinc70"
  egg_minimum_2 = Get time of minimum: point_2, point_3, "Sinc70"
  period = egg_minimum_2 - egg_minimum_1

  if period <= mean_period * 2
    selectObject: "PointProcess degg_smooth"
    degg_maximum_point_1 = Get nearest index: egg_minimum_1
    degg_maximum = Get time from index: degg_maximum_point_1

    if degg_maximum <= egg_minimum_1
      degg_maximum = Get time from index: degg_maximum_point_1 + 1
    endif

    selectObject: "Sound degg_smooth"
    degg_minimum = Get time of minimum: degg_maximum, egg_minimum_2, "Sinc70"

    degg_maximum_rel = (degg_maximum - egg_minimum_1) / period
    degg_minimum_rel = (degg_minimum - egg_minimum_1) / period

    time = egg_minimum_1 - start
    proportion = (egg_minimum_1 - start) / (end - start)

    result_line$ = "'speaker$', 'filename$', 'date$', 'stimulus$', 'egg_minimum_1',
      ... 'time', 'proportion', 'degg_maximum_rel', 'degg_minimum_rel'"

    appendFileLine: "'result_file$', "'result_line$"
  endif
endfor

```

“smoothing”

```

procedure smoothing : .width
  .weight = .width / 2 + 0.5

  .formula$ = "( "

  for .w to .weight - 1
    .formula$ = .formula$ + string$(.w) + " * (self [col - " + string$(.w) + "] +
      ...self [col - " + string$(.w) + "]) + "
  endfor

  .formula$ = .formula$ + string$(.weight) + " * (self [col]) ) / " +
    ...string$(.weight ^ 2)

```

```
Formula: .formula$  
endproc
```

Word DEGG tracing

degg-tracing-word.praat

```
<<<smoothing>>>  
  
<<<get files list word>>>  
  
<<<file loop word>>>
```

“get files list word”

```
form dEGG tracing  
  word project pilot  
  word speaker sc01  
  comment Specify the lower and upper frequency (in Hz) for filtering:  
  real lower 40  
  real upper 10000  
  comment Specify the smooth width "m" (the number of points):  
  real smooth_width 11  
endform  
  
directory$ = "../data/derived/egg/'speaker$'"  
directory_textgrid$ = "../data/derived/ultrasound/'speaker$'/audio"  
  
result_file$ = "../results/'speaker$'-degg-tracing-word.csv"  
header$ = "speaker,file,word,time,rel.time,proportion,maximum,minimum"  
writeFileLine: "'result_file$'", "'header$'"  
  
Create Strings as file list: "filelist", "'directory$'/*.wav"  
files = Get number of strings
```

“file loop word”

```
for file to files  
  selectObject: "Strings filelist"  
  file$ = Get string: file  
  filename$ = file$ - ".wav"  
  
  Read Strings from raw text file: "'directory_textgrid$'/'filename$'.txt"  
  prompt$ = Get string: 1  
  stimulus$ = extractWord$(prompt$, " ")
```

```

    Read separate channels from sound file: "'directory$'/'file$'"

    Read from file: "'directory_textgrid$'/'filename$'.TextGrid"

    start = Get starting point: 2, 2
    end = Get end point: 2, 2

    selectObject: "Sound 'filename$'_ch2"
    Extract part: start, end, "rectangular", 1, "yes"
    Rename: "egg"

    <<<main function>>>
endfor

```

Get durations

get-durations.praat

```

form Get vowel duration
    word project pilot
    word speaker sc01
    comment Supported languages: it, pl
    word language it
endform

if language$ == "it"
    label_lang$ = "dico"
elif language$ == "pl"
    label_lang$ = "mowie"
else
    exit "The language you selected is not valid"
endif

directory$ = "../data/derived/ultrasound/'speaker$'/alignment"

result_file$ = "../results/'speaker$'-vowel-durations.csv"

header$ = "index,speaker,file,word,time,word.duration,c1.duration,vowel.duration,
    ...closure.duration,rvot,c2.duration,v2.duration,sentence.duration"
writeFileLine: result_file$, header$

bursts = Read from file: "'directory$'/'speaker$'-burst.TextGrid"

palgin = Read from file: "'directory$'/'speaker$'-palgin.TextGrid"

intervals = Get number of intervals: 2

fileNames = Read from file: "'directory$'/'speaker$'-filenames.TextGrid"

```



```

index = 0

for interval to intervals
  selectObject: palign
  label$ = Get label of interval: 2, interval
  if label$ == label_lang$
    index += 1
    word$ = Get label of interval: 2, interval + 1
    start_target = Get start time of interval: 2, interval + 1
    end_target = Get end time of interval: 2, interval + 1
    word_duration = (end_target - start_target) * 1000
    start_consonant = Get interval at time: 1, start_target
    start_vowel = Get start time of interval: 1, start_consonant + 1
    c1_duration = (start_vowel - start_target) * 1000
    end_vowel = Get end time of interval: 1, start_consonant + 1
    end_consonant2 = Get end time of interval: 1, start_consonant + 2
    v_duration = (end_vowel - start_vowel) * 1000
    v2_duration = (end_target - end_consonant2) * 1000
    sentence_interval = Get interval at time: 3, start_target
    start_sentence = Get start time of interval: 3, sentence_interval
    end_sentence = Get end time of interval: 3, sentence_interval
    sentence_duration = end_sentence - start_sentence

    selectObject: bursts
    burst_interval = Get nearest index from time: 1, end_vowel
    burst = Get time of point: 1, burst_interval
    if burst < end_vowel or burst > end_sentence
      burst = undefined
    endif

    closure = (burst - end_vowel) * 1000
    rvot = (end_consonant2 - burst) * 1000
    consonant_duration = closure + rvot

    selectObject: fileNames
    fileName = Get interval at time: 1, start_vowel
    fileName$ = Get label of interval: 1, fileName

    result_line$ = "'index','speaker$','fileName$','word$','start_target',
      ...'word_duration','c1_duration','v_duration','closure','rvot',
      ...'consonant_duration','v2_duration','sentence_duration'"
    appendFileLine: "'result_file$'", "'result_line$'"
  endif
endfor

removeObject: palign, bursts

```

Burst detection

This script detects the burst in the consonant following the target vowels (C2). The algorithm is based on @avanthapadmanabha2014.

burst-detection.praat

```
<<<get alignment>>>
```

```
<<<find consonant>>>
```

We start by identifying the interval that corresponds to C2.

“find consonant”

```
speech_intervals = Get number of intervals: 3
sound = Read from file: "'directory_alignment$'/'speaker$'.wav"
textgrid = To TextGrid: "burst","burst"

for speech_interval to speech_intervals
  selectObject: palign
  speech_label$ = Get label of interval: 3, speech_interval
  if speech_label$ == "speech"
    speech_start = Get start time of interval: 3, speech_interval
    token_interval = Get interval at time: 2, speech_start
    token_end = Get end time of interval: 2, token_interval
    phone_interval = Get interval at time: 1, token_end
    start_consonant = Get start time of interval: 1, phone_interval + 2
    end_consonant = Get end time of interval: 1, phone_interval + 2

    selectObject: sound
    sound_consonant = Extract part: start_consonant, end_consonant,
      ..."rectangular", 1, "yes"

    <<<filter>>>

    <<<plosion index>>>

    selectObject: textgrid
    if burst <> undefined
      Insert point: 1, burst, "burst"
    endif
  endif
endfor

selectObject: textgrid
Save as text file: "'directory_alignment$'/'speaker$'-burst.TextGrid"
```

To calculate the plosion index, it is first necessary to filter the sound file.

“filter”

```
Filter (pass Hann band): 400, 0, 100
sound_band = selected("Sound")

spectrum = To Spectrum: "no"
Rename: "original"

spectrum_hilbert = Copy: "hilbert"
Formula: "if row=1 then Spectrum_original[2,col] else -Spectrum_original[1,col] fi"
sound_hilbert = To Sound
samples = Get number of samples
Formula: "abs(self)"
matrix = Down to Matrix
period = Get column distance
```

We can now calculate the plosion index.

“plosion index”

```
m1_time = 0.006
m2_time = 0.016

for sample from 1 to samples
  current = sample * period
  selectObject: sound_hilbert
  mean_before = Get mean: 1, current - m1_time - m2_time, current - m1_time
  mean_after = Get mean: 1, current + m1_time, current + m1_time + m2_time
  window_average = (mean_before + mean_after) / 2
  current_value = Get value at time: 1, current, "Sinc70"
  plosion = current_value / window_average

  if plosion == undefined
    plosion = 0
  elif plosion < 3
    plosion = 0
  endif

  selectObject: matrix
  Set value: 1, sample, plosion
endfor

To Sound
Shift times by: start_consonant
To PointProcess (extrema): 1, "yes", "no", "Sinc70"
half_consonant = start_consonant + ((end_consonant - start_consonant) / 3) * 2
Remove points between: start_consonant, half_consonant
burst = Get time from index: 1
```

Get measurements

This script extracts several durations related to voicing. The main function `merge` is a loop that reads the TextGrids from the derived ultrasound and EGG folders and merges the tier with the gestures from the ultrasound and the tier with the voiced/unvoiced intervals from the EGG.

get-measurements.praat

```
<<<read>>>
```

```
<<<merge>>>
```

This is the form that prompts the user to input the directories of the derived ultrasound (`directory_us`) and EGG (`directory_egg`) data, and the ID of the participant (`speaker`). Do not include the participant folder in the path because it will be automatically included in the main function.

“read”

```
form Get measurements
  word speaker sc01
endform

directory_us_annotations$ = "../data/derived/ultrasound/'speaker$'/
  ...annotations"
directory_egg_vuv$ = "../data/derived/egg/'speaker$'"

createDirectory("../data/derived/merged/'speaker$'")
directory_out$ = "../data/derived/merged/'speaker$'"

result_file$ = "../results/'speaker$'-measurements.csv"
result_header$ = "speaker,word,target,max,release,voff,voffr"
writeFileLine: result_file$, result_header$

Create Strings as file list: "filelist_us", "'directory_us_annotations$'/*.TextGrid"
files_us = Get number of strings

Create Strings as file list: "filelist_egg", "'directory_egg_vuv$'/*.TextGrid"
files_egg = Get number of strings
```

“merge”

```
for file from 1 to files_us
  selectObject: "Strings filelist_us"
  file$ = Get string: file
  Read from file: "'directory_us_annotations$'/'file$'"
```

```

filename$ = selected$("TextGrid")

num_tiers = Get number of tiers

if num_tiers == 4
    Extract one tier: 4

    selectObject: "Strings filelist_egg"
    Read from file: "'directory_egg_vuv$'/'filename$'-vuv.TextGrid"

    selectObject: "TextGrid PointTier_0"
    plusObject: "TextGrid " + filename$ + "-vuv"

    Merge

    Set tier name: 1, "gestures"
    Insert interval tier: 3, "stimulus"

    Read Strings from raw text file: "'directory_us_annotations$'/'filename$'.txt"
    prompt$ = Get string: 1
    stimulus$ = extractWord$(prompt$, " ")

    selectObject: "TextGrid merged"
    Set interval text: 3, 1, stimulus$

    Save as text file: "'directory_out$'/'filename$'-merged.TextGrid"

    <<<calculate>>>
endif
endfor

```

For the current TextGrid, get the number of points in the `gestures` point tier and, if `number_of_points > 0`, loop through the points. If the point is labelled `target_TT` or `target_TD`, get the time and save it to `target`. Else, write an empty value to `target`, and if the label is `max_TT` or `max_TD`, get the time and write it to `max`. Else, write an empty to `max`, and if the label is `release_TT` or `release_TD`, write the value to `release`. Else, write an empty to `release`.

“calculate”

```

number_of_points = Get number of points: 1

target = undefined
max = undefined
release = undefined
voff = undefined
voffr = undefined

if number_of_points > 0
    for point to number_of_points

```

```

point_label$ = Get label of point: 1, point
if point_label$ == "target_TT" or point_label$ == "target_TD"
    target = Get time of point: 1, point
    vuv = Get interval at time: 2, target
    vuv_label$ = Get label of interval: 2, vuv
    if vuv_label$ == "U"
        voff = Get starting point: 2, vuv
    else
        voffr = 0
    endif
elseif point_label$ == "max_TT" or point_label$ == "max_TD"
    max = Get time of point: 1, point
    if target == undefined
        vuv = Get interval at time: 2, max
        vuv_label$ = Get label of interval: 2, vuv
        if vuv_label$ == "U"
            voff = Get starting point: 2, vuv
        else
            voffr = 0
        endif
    endif
elseif point_label$ == "release_TT" or point_label$ == "release_TD"
    release = Get time of point: 1, point
endif
endfor
if voffr <> 0
    if voff == undefined or release == undefined
        voffr = undefined
    else
        voffr = (release - voff) * 1000
    endif
endif
endif

result_line$ = "'speaker$', 'stimulus$', 'target', 'max', 'release',
... 'voff', 'voffr'"
appendFileLine: result_file$, result_line$

```

Get the number of a tier based on the name

The following is a procedure that returns the number of a tier in a TextGrid given the name of that tier. The value is returned to `getTierNumber.return`.

“get tier number”

```

procedure getTierNumber: .tierName$
    .numberOfTiers = Get number of tiers

```

```

        .index = 1
    repeat
        .current$ = Get tier name: .index
        .index += 1
    until .current$ == .tierName$ or .index > .numberOfTiers
    if .index > .numberOfTiers
        exitScript: "The selected TextGrid does not have a tier named '".tierName$'.'"
    else
        .return = .index - 1
    endif
endproc

```

Get duration of voicing in vowels

voicing-duration.praat

```
<<<voicing setup>>>
```

```
<<<voicing loop>>>
```

“voicing setup”

```

form Get duration of voicing
    word speaker sc01
endform

vuvDirectory$ = "../data/derived/egg/'speaker$"
palginDirectory$ = "../data/derived/ultrasound/'speaker'/audio"
resultsFile$ = "../results/'speaker$'-voicing.csv"
resultsHeader$ = "index,speaker,file,rec.date,word,voicing.start,voicing.end,voicing.duration,sente
writeFileLine: resultsFile$, resultsHeader$

Create Strings as file list: "vuvList", "'vuvDirectory$'/*.TextGrid"
numberOfVuv = Get number of strings
index = 0

```

“voicing loop”

```

for vuv to numberOfVuv
    selectObject: "Strings vuvList"
    vuvFile$ = Get string: vuv
    vuvTextGrid = Read from file: "'vuvDirectory$'/'vuvFile$'"
    vuvTextGrid$ = selected$("TextGrid")
    palginTextGrid$ = vuvTextGrid$ - "-vuv"

    Read Strings from raw text file: "'palginDirectory$'/'palginTextGrid$'.txt"

```

```

recDate$ = Get string: 2

palignTextGrid = Read from file: "'palignDirectory$'/'palignTextGrid$'-palign.TextGrid"
plusObject: vuvTextGrid
Merge
numberOfWords = Get number of intervals: 3

<<<words loop>>>
endfor

```

“words loop”

```

for word to numberOfWords
word$ = Get label of interval: 3, word
if word$ == "dico" or word$ == "mowie"
    index = index + 1
    wordStart = Get start time of interval: 3, word + 1
    segment = Get interval at time: 2, wordStart
    vowelStart = Get start time of interval: 2, segment + 1
    vowelEnd = Get end time of interval: 2, segment + 1
    midPoint = vowelStart + (vowelEnd - vowelStart)
    voiced = Get interval at time: 1, midPoint
    voicedStart = Get start time of interval: 1, voiced
    voicedEnd = Get end time of interval: 1, voiced
    voicing = (voicedEnd - voicedStart) * 1000
    stimulus$ = Get label of interval: 3, word + 1

    sentenceInterval = Get interval at time: 4, vowelStart
    sentenceStart = Get start time of interval: 4, sentenceInterval
    sentenceEnd = Get end time of interval: 4, sentenceInterval
    sentenceDuration = sentenceEnd - sentenceStart

    resultLine$ = "'index','speaker$','palignTextGrid$','recDate$','stimulus$',
        ...'voicedStart','voicedEnd','voicing','sentenceDuration'"
    appendFileLine: resultsFile$, resultLine$
endif
endfor

```

Headers

“sync header”

```

#####
# sync-egg.praat v1.0.0
#####
# MIT License
#

```



```

# Copyright (c) 2016 Stefano Coretta
#
# Permission is hereby granted, free of charge, to any person obtaining a copy
# of this software and associated documentation files (the "Software"), to deal
# in the Software without restriction, including without limitation the rights
# to use, copy, modify, merge, publish, distribute, sublicense, and/or sell
# copies of the Software, and to permit persons to whom the Software is
# furnished to do so, subject to the following conditions:
#
# The above copyright notice and this permission notice shall be included in all
# copies or substantial portions of the Software.
#
# THE SOFTWARE IS PROVIDED "AS IS", WITHOUT WARRANTY OF ANY KIND, EXPRESS OR
# IMPLIED, INCLUDING BUT NOT LIMITED TO THE WARRANTIES OF MERCHANTABILITY,
# FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT. IN NO EVENT SHALL THE
# AUTHORS OR COPYRIGHT HOLDERS BE LIABLE FOR ANY CLAIM, DAMAGES OR OTHER
# LIABILITY, WHETHER IN AN ACTION OF CONTRACT, TORT OR OTHERWISE, ARISING FROM,
# OUT OF OR IN CONNECTION WITH THE SOFTWARE OR THE USE OR OTHER DEALINGS IN THE
# SOFTWARE.
#####
# This script syncs the audio files acquired by the Laryngograph with the audio
# files exported from AAA. Syncing is obtained through pair-wise
# cross-correlation of the audio files. The cross-correlation function returns
# the off-set in seconds between two files. The off-set is used to remove the
# leading audio from the longer file.
#
# Input: - .wav stereo files from the Laryngograph recordings (ch1 = audio, ch2 =
# EGG), saved in a folder
#       - .wav mono files exported from AAA, saved in a separate folder
# Output: - .wav stereo files (ch1 = audio, ch2 = EGG) whose start time is
# synced with the start time of the correspondet AAA file
#####

```